



TRITTOOLS CO

Where Discovery Knows No Bounds

Precision Tweezers





Page 3	PT-00 - PT-00C - PT-0 - PT-0C9 - PT-0C10 - PT-0C11
Page 4	PT-0A - PT-0R - PT-1 - PT-2 - PT-2A - PT-2AB
Page 5	PT-3 - PT-3C - PT-4 - PT-4A - PT-5 - PT-5A
Page 6	PT-51S - PT-6 - PT-7 - PT-7A - PT-7B - PT-7E
Page 7	PT-8 - PT-26 - PT-27 - PT-S - PT-SA - PT-SB
Page 8	PT-SS - PT-P - PT-F - PT-H - PT-PC - PT-107
Page 9	PT-114 - PT-231 - PT-232 - PT-233 - PT-35A - PT-36A
Page 10	PT-37S - PT-29D - PT-30D - PT-190M - PT-A47 - PT-A48
Page 11	PT-A49 - PT-29 - PT-15A - PT-15AP - PT-571 - PT-572
Page 12	PT-573 - PT-574 - PT-53D - PT-4D - PT-5D - PT-6D
Page 13	PT-4DW - PT-5DW - PT-6DWF - PT-AA - PT-MM - PT-GG
Page 14	PT-HH - PT-NN - PT-OO - PT-PP - PT-R - PT-RR
Page 15	PT-W - PT-DD - PT-J05 - PT-J06 - PT-J07 - PT-J08
Page 16	PT-32 - PT-G01 - PT-G02 - PT-G03 - PT-G04 - PT-G05
Page 17	PT-G06 - PT-G07 - PT-001 - PT-002 - PT-C1 - PT-C2
Page 18	PT-C3 - PT-C4 - PT-C5 - PT-C6 - PT-C7 - PT-C8
Page 19	PT-C9 - PT-C10 - PT-C11 - PT-C12 - PT-C13 - PT-C14
Page 20	PT-C15





Precision Tweezers



Thick and Flat, very strong points

Size	Art No
------	--------

120mm	PT-00
-------	-------



Thick and Flat, very strong points

Size	Art No
------	--------

100mm	PT-00C
-------	--------



Flat, fine points

Size	Art No
------	--------

120mm	PT-0
-------	------



Flat, fine points

Size	Art No
------	--------

90mm	PT-0C9
------	--------



Flat, fine points

Size	Art No
------	--------

100mm	PT-0C10
-------	---------




Flat, fine points

Size	Art No
------	--------

110mm	PT-0C11
-------	---------



Precision Tweezers



Fine, for watch repairers

Size

Art No

120mm

PT-0A




Pointes fines, action reversible

Size

Art No

120mm

PT-0R



Heavy, sharp

Size

Art No

120mm

PT-1



Strong, fine points

Size

Art No

120mm

PT-2



Flat, rounded points, flache backen

Size

Art No

120mm

PT-2A



Tips as 2A but curved

Size

Art No

120mm

PT-2AB



Precision Tweezers



Very sharp

Size

Art No

120mm

PT-3



Short, very sharp

Size

Art No

110mm

PT-3C



Extra fine points

Size

Art No

110mm

PT-4



Stronger points

Size

Art No

110mm

PT-4A



Super fine points

Size

Art No

110mm

PT-5



Oblique points

Size

Art No

115mm

PT-5A



Precision Tweezers



Super fine points

Size

Art No

115mm

PT-51S



Sharp hooked points

Size

Art No

115mm

PT-6



Fine curved points

Size

Art No

110mm

PT-7



Strong curved points

Size

Art No

110mm

PT-7A



Serrated points

Size

Art No

110mm

PT-7B



Flat, fine points

Size

Art No

110mm

PT-7E



Precision Tweezers



For straightening pivots and wires

Size

Art No

115mm

PT-8



Very strong, sharp points

Size

Art No

135mm

PT-26



Strong, fine points

Size

Art No

135mm

PT-27



Plain, fine points

Size

Art No

120mm

PT-S



For everyday use, strong points

Size

Art No

120mm

PT-SA



For everyday use, fine points

Size

Art No

120mm

PT-SB



Precision Tweezers



Very long and narrow, fine points

Size

Art No

135mm

PT-SS



Plain, strong points

Size

Art No

120mm

PT-P



For flattening hands, flat end

Size

Art No

120mm

PT-F



Wide shanks, fine

Size

Art No

90mm

PT-H



Tweezers for laying contacts

Size

Art No

90mm

PT-PC



Strong and serated points

Size

Art No

107mm

PT-107



Precision Tweezers



Strong and tapering points

Size

Art No

114mm

PT-114



Cheap tweezers, serrated points

Size

Art No

120mm

PT-231



Heavy gripping tweezers

Size

Art No

160mm

PT-232



Heavy gripping tweezers, slightly angled jaws

Size

Art No

160mm

PT-233



Tweezers for flat components, wafers and disc shaped condensers

Size

Art No

120mm

PT-35A



As 35A but with angled points

Size

Art No

120mm

PT-36A



Precision Tweezers



Wafer handling tweezers with stop pins

Size

Art No

120mm

PT-37S



Very useful when removing dipoles

Size

Art No

152mm

PT-29D



Very useful when removing dipoles but with bent points

Size

Art No

152mm

PT-30D



Fiber glass tweezers, acid resistant

Size

Art No

110mm

PT-190M



Carbon coated tweezers carbon = red

Size

Art No

PT-A47

Inbox = green

Size

Art No

PT-A48



Precision Tweezers



Amagnetic = blue

Size

Art No

PT-A49



For assemblage, with serrated points

Size

Art No

160mm

PT-19



For cutting hairsprings and electronic component wires

Size

Art No

115mm

PT-15A



Same as 15A with parallel cut

Size

Art No

115mm

PT-15AP



Fit components of .098.196 inch for positioning and desoldering

Size

Art No

145mm

PT-571



Fit components of .098.196 inch for positioning and desoldering

Size

Art No

145mm

PT-572



Precision Tweezers



Fit components of .098.196 inch for positioning and desoldering

Size	Art No
------	--------

145mm	PT-573
-------	--------



Tweezer for assembling and welding flat electronic parts

Size	Art No
------	--------

145mm	PT-574
-------	--------



Wafer tweezers with flat, lower paddle 3 teeth

Size	Art No
------	--------

130mm	PT-3D
-------	-------



Same shape as 3D with 4 teeth

Size	Art No
------	--------

135mm	PT-4D
-------	-------



Same shape as 3D with 5 teeth

Size	Art No
------	--------

130mm	PT-5D
-------	-------



Same shape as 3D with 6 teeth

Size	Art No
------	--------

130mm	PT-6D
-------	-------



Precision Tweezers



Wafer tweezer, no teeth

Size

Art No

130mm

PT-4DW



Wafer tweezer, no teeth

Size

Art No

130mm

PT-5DW



Wafer tweezer with graduated lower paddle and 6 teeth

Size

Art No

130mm

PT-6DWF



Strong tweezers

Size

Art No

125mm

PT-AA



Fine beveled points

Size

Art No

125mm

PT-MM



Large and fine points

Size

Art No

130mm

PT-GG



Precision Tweezers



Very strong, for ordinary use

Size

Art No

115mm

PT-HH



Tweezer for semi-conductors

Size

Art No

115mm

PT-NN



Light and slender, for ordinary use

Size

Art No

110mm

PT-00



Tweezers for general

Size

Art No

115mm

PT-PP



Long Blades, tapering points

Size

Art No

140mm

PT-R



Medium Points long and strong blades

Size

Art No

140mm

PT-RR



Precision Tweezers



Extra heavy wide blades, medium fine points

Size

Art No

160mm

PT-W



Heavy and large blades, medium points

Size

Art No

160mm

PT-DD



Soldering tweezers straight, reverse action

Size

Art No

160mm

PT-J05



Soldering tweezers straight with insulated handle

Size

Art No

165mm

PT-J06



Soldering tweezers straight with insulated handle

Size

Art No

165mm

PT-J07



Soldering tweezers straight with insulated handle

Size

Art No

170mm

PT-J08



Precision Tweezers



Soldering tweezers

Size

Art No

200mm

PT-32



Jewellers tweezers with horizontal lines

Size

Art No

140mm

PT-G01



Jewellers tweezers with horizontal lines

Size

Art No

160mm

PT-G02



Jewellers tweezers with horizontal lines

Size

Art No

165mm

PT-G03



Jewellers tweezers with locking pin

Size

Art No

140mm

PT-G04



Jewellers tweezers with locking pin

Size

Art No

150mm

PT-G05



Precision Tweezers



Jewellers tweezers with vertical lines

Size

Art No

160mm

PT-G06



Gams tweezers with shovel

Size

Art No

160mm

PT-G07



Universal tweezers with guide pin
straight points

Size

Art No

160mm

PT-001



Tweezers with guide pin with curved points

Size

Art No

160mm

PT-002



Tweezer with guide pin with straight points

Size

Art No

150mm

PT-C1



College tweezer with guide pin with curved points

Size

Art No

150mm

PT-C2



Precision Tweezers



Meriam tweezer with guide pin and double curved points

Size

Art No

160mm

PT-C3



General purpose thumb tweezer

Size

Art No

160mm

PT-C4



General purpose splinter tweezer

Size

Art No

160mm

PT-C5



College tweezer coated with guide pin with curved points

Size

Art No

160mm

PT-C6



Meriam tweezer coated with guide pin and double curved points

Size

Art No

150mm

PT-C7



Tweezer with extra opening guide pin and curved points

Size

Art No

150mm

PT-C8



Precision Tweezers



Titanium tweezers very sharp short

Size

Art No

120mm

PT-C9



Titanium tweezers very sharp

Size

Art No

120mm

PT-C10



Titanium tweezers extra fine points

Size

Art No

120mm

PT-C11



Titanium tweezers super fine points

Size

Art No

110mm

PT-C12



Titanium tweezers oblique points

Size

Art No

115mm

PT-C13



Titanium tweezers flat fine points

Size

Art No

110mm

PT-C14



Precision Tweezers



Tweezer coated very long and narrow fine points

Size

Art No

135mm

PT-C15



Plastic Tweezer Serrated with English Pattern Printed Handle

Size

Art No

12.5cm

PT-C21




Plastic Tweezer Serrated with Printed Handle

Size

Art No


11 cm
13.5cm

PT-C24
PT-C25




Size

Art No



Size

Art No



Size

Art No



TRITOOOLS CO

Where Discovery Knows No Bounds

.Contact

Factory Address

Near Round About S.I.E, Sialkot-Pakistan

Mobile/WhatsApp

0092 307 9695921

Email

export@tritoolsco.com

Website

www.tritoolsco.com

